

Farmasötik-Hızlı Test

Uzunluk (mm)	Kapasite (ml)	Model	Numarası	GTIN
75	0.2	Farmasötik- Hızlı Test	120318	6260807503669
140	1.5	Dereceli	120328	6260807503546
156	2	Dereceli	120329	6260807503553
161	3	Dereceli	120330	6260807503560
260	7.5	Ekstra uzun	120325	6260807503638
280	7.5	Ekstra uzun	120326	6260807503645
300	7.5	Ekstra uzun	120327	6260807503652
100	4	Körüklü	120319	6260807503577
120	4	Körüklü	120320	6260807503584
140	4	Körüklü	120321	6260807503591
176	5.5	Körüklü	120322	6260807503607
196	5.5	Körüklü	120323	6260807503614
216	5.5	Körüklü	120324	6260807503621

+98 21 88545922-9

+98 21 88765561

+98 9 123340197

@Poleideal

Poleideal

www.medpip.com



Pole Ideal Pars

Ostad Motahari Caddesi, Mofatteh Caddesi

sonrası, Jahantab Caddesi, Naghdi

Caddesi, No:12, Kat:1,Tehran,İran

Posta kodu: 1576635714 Posta kutusu: 15875-9483

Çat.No:00329



Pasteur Pipette

Pastör pipetleri
Pole Ideal Pars





İhtiyaçlarınıza göre seçebilirsiniz.



Kapasite (ml)	0.2	1.5	2	3	4	4	4	5.5	5.5	5.5	7.5	7.5	7.5
Uzunluk (mm)	75	140	154	161	100	120	140	176	196	216	260	280	300
φ Ampul (mm)	8	13	11	13	21	21	21	24	24	24	19	19	19
Numarası	120318	120328	120329	120330	120319	120320	120321	120322	120323	120324	120325	120326	120327



www.medpip.com



Pole Ideal Pars